## SUPRABONE + ASC Kit: KEMİK İLİĞİ ASPIRASYON VE UYGULAMA KİTİ İLE KEMİK GREFT MATERYALİ KULLANIM TALİMATI

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 **bmtcalsis**  
implanting new ideas

CE  
1984

# SUPRABONE + ASC KIT: BONE GRAFT MATERIAL AND BONE MARROW ASPIRATION AND APPLICATION KIT INSTRUCTIONS FOR USE

## Classification

Suprabone + ASC Kit is a medical device, not a drug.

## Material

Suprabone + ASC Kit consists of-TCP granules, a needle and a syringe. Suprabone - TCP is pure  $\beta$ -tricalcium phosphate ( $\beta$ -TCP). When Suprabone - TCP is implanted it acts as a biocompatible and resorbable osteoconductive support matrix and supports new bone formation at the defect site. Over time, Suprabone - TCP is resorbed and replaced by newly formed bone. As it does not contain any human or animal origin materials it does not pose a threat of disease transmission. The product contains silicate bioglass additive for enhanced bone bonding. As for the ASC Kit, it is composed of a sterile needle and syringe which is used to aspirate bone marrow.

## Indications-For-Use

The implantable component of Suprabone + ASC Kit, Suprabone - TCP, is an osteoconductive bone graft material in the form of granules. It is indicated for defects in non-load bearing regions such as the limbs, spine, pelvis and cranium. The granules and powder can be molded to fit the defect shape and size and then pressed in to the defect site. The bone marrow aspirate obtained with the ASC Kit may be mixed with the Suprabone - TCP granules or with other bone graft substitute materials. Suprabone - TCP functions as a temporary tissue scaffold and is not intended to provide mechanical support throughout the treatment period. The implant is biocompatible and radiopaque and will be resorbed in a controlled manner within 24 months.

## Contraindications

Suprabone + ASC Kit is not indicated for any other uses other than those stated here. It should not be used in the presence of any contraindication/s. Suprabone + ASC Kit should not be used to provide structural support in the skeletal system or to gain screw fixation.

Suprabone + ASC Kit must not be used in patients with any of the following conditions:

- Existing acute or chronic infections, particularly at the implantation site
- Severe vascular or neurological disease
- Poorly vascularized implantation site
- Uncontrolled diabetes
- Severe degenerative diseases
- Hypercalcemia, abnormal calcium metabolism

- Inflammatory bone disease
- Dispersed malignant tumors
- Severe renal dysfunction
- Open epiphyseal plated (pediatric patients)
- Who are uncooperative and who cannot or will not post-operative instructions, including those who display drug and/ or alcohol abuse

## Precautions

Suprabone + ASC Kit must be used by surgeons familiar with bone repair and replacement techniques. **SUPRABONE + ASC KIT IS NOT INTENDED FOR LOAD BEARING APPLICATIONS.** To provide structural support and to keep the implant in an environment free from loading, it is important to mechanically secure the area around the implantation site with rigid fixation. To achieve the best bone regeneration results there should be **MAXIMUM SURFACE AREA CONTACT** between the implant and the surrounding bone tissue. Suprabone + ASC Kit must not be used to gain screw fixation. The effect of Suprabone + ASC Kit on patients with the following conditions is unknown:

- Long term infection
- Pregnancy and nursing
- Bone radiation therapy
- Cardiovascular disease
- Metabolic bone disease
- Renal disease

The effects of Suprabone + ASC Kit bone graft material on pediatrics patients and combinations with other materials (e.g. antibiotics or serum) is unknown.

## Possible Complications

Reoperation may be required to remove or replace the implant in specific medical cases or when device failure occurs. Possible complications include but are not limited to:

- Wound complications including local immunological reactions, hematoma, swelling and fluid accumulation, edema, tissue thinning, infection, bone fracture, and other complications that are possible with any surgery
- Disintegration of the implant with or without generation of particulate debris due to the application of a load
- Bone deformity and loss of contour at the defect site
- Allergic reaction to the product

## Warnings

Suprabone -TCP is sold inside a double-blistered closed packet. The ASC Kit comes inside a separate packet within the product box. The product is sterile. **DO NOT USE IF THE PACKAGE INTEGRITY OF ANY COMPONENT HAS BEEN COMPROMISED.**

Read expiry date before use and **DO NOT USE IF THE EXPIRATION DATE HAS BEEN EXCEEDED.** The inner blister should be opened immediately prior to implantation in a sterile environment.

Since Suprabone appears moderately opaque in

X-rays, areas surrounding the implant may be masked on radiographic images. A suitable method of immobilization of the immobilization of Suprabone material in the defect area is required to prevent potential migration and the material should only be used in surgical procedures where bone grafts are adequately supplied.

Since both hard tissues and Suprabone have very short T2 values, it may be difficult to visualize the product with conventional MRI.

Suprabone + ASC Kit is single-use only. **DO NOT RE-STERILIZE AND RE-USE THE PRODUCT AND THE ASC KIT.**

## Application

**Step1:** Open the outer blister in the operating theatre and transfer the inner blister to a sterile area. Open the inner blister in a sterile area prior to implantation. **Step2:** Implant the material. The material can be molded to fit the defect site and mixed with the bone marrow aspirate. The material can be gently and carefully pressed into the defect site. **SUPRABONE - TCP SHOULD BE IN DIRECT CONTACT WITH ALL SURFACES WITHIN THE DEFECT.**

**Step3:** To prevent the implant from moving or migrating secure the surgical site after implantation. If there is an excess of fluid in the surgical area appropriate actions can be taken to stop the bleeding (e.g. cauterization, aspiration or application of bone wax). If the material is not implanted successfully remove the implant and start again with a new dose of Suprabone - TCP

## Application of the ASC (Autologous Stem Cell) Kit

An appropriate amount of bone marrow aspirate is taken from the iliac crest or vertebral body to be mixed with Suprabone - TCP and applied as mentioned above.

## Storage

Store in a **DRY PLACE AT ROOM TEMPERATURE.** Optimal storage conditions: **+5°C - +39°C (41-102.2°F)**, <70% relative humidity. Direct contact with heating systems or storage under direct sunlight should be avoided.

## Shelf Life and Disposal

The expiration date is printed on the label. **DO NOT USE AFTER EXPIRATION DATE.** Suprabone + ASC Kit is environment - friendly. No special disposal procedure is necessary. Packaging material is recyclable oid direct contact with heating systems or storage under direct sunlight.

## Ordering Information

Suprabone + ASC Kit is a sterile osteoconductive bone graft material. Suprabone + ASC Kit is presented with a detailed instruction for use. Suprabone - TCP granules have been sterilized with gamma radiation and the ASC Kit has been sterilized with ethylene oxide treatment. Suprabone + ASC Kit bone filler material, along with its mixing container, is provided inside a double blister packet for transport and storage. In addition to this booklet, there are labels inside the box for patient documentation. For further information on the product and its uses, please contact **BMT Calsis Sağlık Teknolojileri San. Tic. A.Ş.**

## Symbol Key for Product Box and Label

 Manufacturer

 Date of Manufacture

 Use-by Date

 Batch Code

 Catalogue Number

 Sterilized Using Irridation

 Sterilized using ethylene oxide

 Do not Resterilize

 Do not Use if Package is Damaged

 Do not re-use

 Keep away from sunlight

 Keep dry

 <sup>39°C</sup> Temperature limit

 Consult instructions for use

 CE symbol and Notified Body number

**Please read before use.**

## Manufacturer

BMT Calsis Sağlık Teknolojileri San. Tic. A.Ş.  
Saray Mahallesi 1500 Cad. No:30 A 06980  
Kahramankazan Ankara Türkiye  
Tel: +90 312 405 80 69 Fax: +90 312 440 30 48  
www.bmtcalsis.com

CAL.IFU.03.EN  
Issue Date: 15.03.2019  
Rev01/Rev. Date: 05.04.2019

 **bmtcalsis**  
implanting new ideas

CE  
1984