## SUPRAFELT KIKIRDAK ONARIM MATRİKSİ KULLANIM TALİMATI

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**bmtcalsis**  
implanting new ideas

CE  
1984

## SUPRAFELT CARTILAGE REPAIR MATRIX INSTRUCTIONS FOR USE

### Classification

SupraFelt Cartilage Repair Matrix (CRM) is a Class III medical device, not a drug.

### Material

SupraFelt CRM is a cell free implant material consisting of non-woven pure polyglycolic acid enriched with non-animal originated hyaluronic acid. SupraFelt is a cartilage repair matrix that is biocompatible and bioresorbable and supports new cartilage formation in the defect area when implanted. SupraFelt resorbs while it is replaced by the new cartilage tissue formed. SupraFelt does NOT contain animal or human-derived tissue-derived structures and therefore does NOT carry a risk of transmitting disease. SupraFelt can be implanted by arthroscopic or mini-arthotomy (minimally invasive) methods in local articular cartilage defects, using biodegradable-biocomposite pins or various tissue adhesives. SupraFelt can be implanted by arthroscopic or mini-arthotomy (minimally invasive) methods in local articular cartilage defects, using biodegradable-biocomposite pins or various tissue adhesives. Prior to the implantation, SupraFelt can be elasticized by placing in platelet-rich plasma (PRP), patient's own blood serum or sterile water for at least 2 minutes immediately.

### Indications

SupraFelt CRM is a biodegradable polymeric matrix that repairs articular cartilage. It is indicated to the repair of the articular cartilage defect only. SupraFelt can be implanted into the defect area by arthrotomic or arthroscopic methods also can be shaped in order to fit into the defect area. SupraFelt will act as a temporary scaffold and is not intended to provide mechanical support during the healing process.

### Contraindications

SupraFelt CRM is not indicated for any other use than the ones stated above. SupraFelt should not be used when there is a contraindication. SupraFelt is contraindicated where the device is intended to provide structural support in the skeletal system. Other contraindications are:

- Existing acute or chronic infections, particularly at the implantation site
- Chemotherapy or radiotherapy patients

### Precautions

The effect of SupraFelt on patients with the following conditions is unknown:

- Long term infection
- Pregnancy and nursing
- Radiation bone therapy
- Cardiovascular diseases
- Metabolic bone disease
- Renal diseases

The effects of SupraFelt CRM on pediatric patients and the effects of mixing with other substances (e.g. antibiotics) are also unknown.

### Possible Complications

Re-operation to remove or replace an implant may be required, due to specific medical conditions or device failure.

Possible adverse effects may include, and not limited to:

- Wound complications including hematoma, swelling and fluid accumulation, edema, tissue thinning, infection, bone fracture and other complication that are possible with any surgery.
- Disintegration of the implant with or without generation of particulate debris due to load being applied.
- Cartilage deformation and loss of contour at the site.
- Allergic reaction to the product.

### Warnings

SupraFelt CRM comes in a double blister packages; an inner blister and a sealed aluminum outer pack. SupraFelt is sterile.

### DO NOT USE IF THE PACKAGE INTEGRITY IS LOST.

Read the expiry date before use and **DO NOT USE IF THE EXPIRATION DATE IS EXCEEDED**. Inner blister must be open in a sterile environment just before the implantation.

SupraFelt is a single-use product, **DO NOT RESTERILIZE AND REUSE**. Reusing or resterilisation may cause the fail of the sterility, functionality and shape integrity of the product and may accelerate the material degeneration.

### Application

SupraFelt CRM is only intended for the use of surgeons familiar with and skilled in techniques of cartilage repair or replacement.

**Step1:** Tear the outer aluminum pack and inner blister in a sterile environment just before the implantation.

**Step2:** Coat the SupraFelt by using PRP, patients own serum or sterile water.

**Step3:** After the cartilage defect area is debrided and cleaned microfracture is applied. SupraFelt is implanted with biodegradable/biocomposite stabilizers or various tissue adhesives to cover the entire defect area. The SupraFelt can be shaped

regarding to the shape and size of the defect area.

### Storage

SupraFelt CRM must be kept in its own package between 2 – 25 °C in dry conditions. Direct contact with heating systems or storage under direct sunlight should be avoided.

### Shelf Life and Disposal

The expiration date is printed on the label. **DO NOT USE AFTER THE EXPIRATION DATE**. SupraFelt is environment-friendly. No special disposal procedure needed. Packaging material is recyclable.

### Ordering Information

SupraFelt CRM is a sterile cartilage repair and regenerative matrix. SupraFelt CRM is provided with its detailed instructions for use. This product is completely sterilized via gamma sterilization. SupraFelt CRM is packed in an aluminium pack and it comes with an additional box for transportation and storage. In addition to this instruction for use, there are labels in the box for the patient documentation. For further detailed information about product and instructions, please contact with **BMT Calsis Sağlık Teknolojileri San. Tic. A.Ş.**

### Symbol Key for Product Box and Label

-  Manufacturer
-  Date of Manufacture
-  Use-by Date
-  Batch Code
-  Catalogue Number
-  Sterilized Using Irridation
-  Do not Resterilize
-  Do not Use if Package is Damaged
-  Do not re-use
-  Keep away from sunlight
-  Keep dry
-  Temperature limit
-  Consult instructions for use
-  CE symbol and Notified Body number

**Please read before use.**

  
**Manufacturer**  
BMT Calsis Sağlık Teknolojileri San. Tic. A.Ş.  
Saray Mahallesi 1500 Cad. No:30 A 06980  
Kahramankazan Ankara Türkiye  
Tel: +90 312 405 80 69 Fax: +90 312 440 30 48  
www.bmtcalsis.com

CAL.IFU.04.EN  
Issue Date: 01.02.2019  
Rev01/Rev. Date: 05.04.2019

 **bmtcalsis**  
implanting new ideas

CE  
1984