SUPRABONE - ENJEKTABL PUTTY KEMİK GREFT MATERYALİ KULLANIM TALİMATI

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SUPRABONE PUTTY - INJECTABLE BONE GRAFT SUBSTITUTE INSTRUCTIONS FOR USE

Classification

Suprabone Putty is a medical device (Class III), not a drug.

Material

Suprabone Putty composed of two components [Component A and Component B] for reconstitution prior to application. Component A is a β -tricalcium phosphate (β -TCP) and Component B is a cellulosic carrier. The product contains bioglass additive for enhanced bone bonding.

Suprabone Putty can be prepared in two different forms by combining both components in specified ratios (See Application section below). This feature enables surgeons to prepare injectable or malleable putty form for intended use. It is an injectable osteoconductive support matrix which is biocompatible and resorbable when implanted and it supports new bone formation in the defect area. Suprabone Putty resorbs over time and ultimately replaces with natural bone. Cellulosic carrier (Component B) used in the reconstitution of Suprabone Putty is also biocompatible and degradable and it helps vascularization of the defect area by leaving voids following rapid resorption. Suprabone Putty contains no tissue of human or animal origin therefore carries no risk of disease transmission. Suprabone Putty is a non-self setting material. Suprabone Putty components [Component A and Component B] are sterilized via Gamma Irradiation and Steam, respectively and can be easily applied to the defect area in different forms, i.e. injectable or putty forms (See Application section below).

Indications-For-Use

Suprabone Putty is an injectable fully synthetic, osteoconductive bone substitute material. When it is prepared in another ratio specified given in Application section, it can also be prepared as a malleable putty form. In both forms, it is intended for use in non-load bearing defects in periodontal conditions, extremities, spine, pelvis, cranium, etc. as a filler. Main application areas are cystic cavities, periodontal defects, alveolar fillings, sinus lifting, posterior fusion and anterior fusion with cages in spine surgery, etc. Suprabone Putty may also be combined with other bone grafts (autografts/ allografts) prior to application. Suprabone Putty acts as a temporary scaffold, and is not intended to provide mechanical support during the healing process. The implant is biocompatible, radio-opaque and will be resorbed in a controlled way.

Contraindications

Suprabone Putty is not indicated for any other uses than the ones stated here. It should not be used when there is any contraindication. Suprabone Putty is contraindicated where the device is intended to provide structural support in the skeletal system, and must not be used to gain screw fixation.

Other contraindications are:

- Existing acute or chronic infections, particularly at the implantation site
- Severe vascular or neurological disease
- Poorly vascularized implantation site
- Uncontrolled diabetes
- Severe degenerative diseases
- Hypercalcemia, abnormal calcium metabolism
- Inflammatory bone disease
- Malignant tumors
- Immune deficiency
- Osteomyelitis
- Severely impaired renal function
- Open epiphyseal plates in pediatric patients
- Uncooperative patients who cannot or will not follow post-operative instruction, including individuals who abuse drugs and/or alcohol.

Precautions

Suprabone Putty is only intended for the use of surgeons familiar with, and skilled in techniques of bone repair and replacement. **SUPRABONE PUTTY IS NOT INTENDED FOR LOAD BEARING APPLICATIONS.** It is important to **MAXIMIZE THE CONTACT SURFACE** between existing bone and the implant to ensure proper bone regeneration. Suprabone Putty must not be used to gain screw fixation **DO NOT OVERFILL THE DEFECT AREA AND AVOID LEAKAGE OUT OF THE DEFECT AREA.** Suprabone Putty should be reconstituted prior to use by combining two components.

The effect of Suprabone Putty on patients with the following conditions is unknown:

- Long term infection
- Pregnancy and nursing
- Radiation bone therapy
- Cardiovascular disease
- Metabolic bone disease
- Documented renal disease

The effects of Suprabone Putty in pediatric patients and the effects of mixing with other substances (e.g. antibiotics) are also unknown.

Possible Complications

Reoperation to remove or replace an implant may be required, due to specific medical conditions or device failure. Possible adverse effects may include, and are not limited to:

- Wound complications including hematoma, swelling and fluid accumulation, edema, tissue thinning, infection, bone fracture, and other complications that are possible with any surgery,
- Disintegration of the implant with or without generation of particulate debris due to a load being applied,
- Bone deformity and loss of contour at the site,
- Allergic reaction to the product.

Warnings

Suprabone Putty comes as two components for reconstitution during surgery. Component A and Component B are sterilized in double blister packaging and syringes by using different routes of sterilization. Both components are sterile. **DO NOT USE IF THE PACKAGE INTEGRITY OF ANY COMPONENT (A or B) IS LOST.**


Read expiry date before use and **DO NOT USE IF THE EXPIRATION DATE HAS BEEN EXCEEDED.**

Since Suprabone Putty is opaque to X-rays, areas under or above the implant may be masked on the radiography.

Suitable method of immobilization of Suprabone Putty material in the defect area is required to prevent potential migration, and the material should only be used in surgical procedures where bone grafts are adequately contained. Suprabone Putty is for single use only. **DO NOT ATTEMPT TO RE-STERILIZE OR RE-USE BOTH COMPONENTS.**

Application

Suprabone Putty has two components, namely Component A and Component B. Please follow the instruction below given for Injectable and Putty forms:

Component A

Component B

Step 1: Pull the Tyvek lids off of the outer and inner blisters to expose Component A for reconstitution.

PLEASE RECONSTITUTE THE COMPONENTS IN THE INNER BLISTER USED FOR COMPONENT A.

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Step 2: Tear the outer aluminum pouch and the inner Tyvek pouch of the Component B and inject the gel into the inner blister and mix the ceramic powder (Component A) and gel (Component B) as follow:

2.1 FOR MALLEABLE PUTTY FORM:

2.1.a Use HALF of Component B, mix both components until it becomes homogenous mixture

2.1.b Shape putty material by hand.


2.2 FOR INJECTABLE FORM

2.2.a Use ALL of Component B, mix both components until it becomes homogenous mixture.

2.2.b. Fill the same syringe with the final reconstituted composition.

Step 3: Apply Suprabone Putty into the defect area immediately. **PLEASE DO NOT STORE THE MIXED PRODUCT.** Make sure that **SUPRABONE PUTTY IS IN DIRECT CONTACT WITH ALL SURFACES OF THE DEFECT AREA.**

Step 4: Secure the surgical site after implanting in order to prevent motion, leakage or implant migration. When excess fluid is present in the surgical field, the surgeon may use proper measures (i.e. catheterization, suction and application of bone wax) to reduce bleeding.

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Storage

Store in a **DRY PLACE AT ROOM TEMPERATURE.** Optimal storage conditions: **+5°C - +39°C (41-102,2°F)** <70 relative humidity. Direct contact with heating systems or storage under direct sunlight should be avoided.

Shelf Life and Disposal

The expiration date is printed on the label. **DO NOT USE AFTER THE EXPIRATION DATE.** Suprabone Putty is environment-friendly. No special disposal is necessary.

Ordering Information

Suprabone Putty is a sterile, injectable osteoconductive bone graft substitute reconstituted during surgery. It is provided with detailed instructions-for-use. It has two components (A and B). Both components are sterile (via gamma irradiation and steam). Components of Suprabone Putty bone graft substitute are packed in two different packs. Component A is filled in PET-G blister with Tyvek lids while Component B is provided within PC syringes and packed in double pouches. Both components come within an additional box for transportation and storage. In addition to this booklet, there are labels inside the box for documentation. For further information on the product and its uses, please contact **BMT Calsis Sağlık Teknolojileri San. Tic. A.Ş.**

Symbol Key for Product Box and Label

Manufacturer	Do not Use if Package is Damaged
Date of Manufacture	Do not re-use
Use-by Date	Keep away from sunlight
LOT Batch Code	Keep dry
REF Catalogue Number	5°C - 39°C Temperature limit
STERILE Sterilized Using Irradiation	Consult instructions for use
STERILE Sterilization using steam or dry heat	CE CE symbol and Notified Body number
Do not Resterilize	

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Manufacturer
BMT Calsis Sağlık Teknolojileri San. Tic. A.Ş.
Saray Mahallesi 1500 Cad. No:30 A 06980
Kahramankazan Ankara Türkiye
Tel: +90 312 405 80 69 Fax: +90 312 440 30 48
www.bmtcalsis.com

CAL.IFU.02.EN
Issue Date: 15.03.2019
Rev01/Rev. Date: 05.04.2019

bmtcalsis
implanting new ideas