

Protective Gown Blue



Product Code: **PPE102**



Description

40 gsm PP / SMS fabric - Sewn from light , soft , water resistant two layer nonwoven fabric Elastic knit cuffs

Fastened at the back of the neck with velcro strap

Can decrease the potential for cross -contamination in pandemic response activities when used with mask and gloves

The length is 125 cm but can be produced in different height measurement

Available in Blue colour

Certifications

ISO 13485:2016 / ISO 9001:2015(QMS)
/ CE ATTESTATION OF CONFORMITY
RELATED TO MEDICAL DEVICES
93/42/EEC/CLASS 1