

Protective Gown White



Product Code: **PPE101**

Description

57 gsm NW textile backsheet
- Sewn from light , soft , water
resistant two layer nonwoven
fabric Elastic knit cuffs

Fastened at the back of the neck
with velcro strap

Can decrease the potential
for cross -contamination in
pandemic response activities
when used with mask and gloves

The length is 125 cm but can
be produced in different height
measurement

Available in White colour

Certifications

EN 13795-1:2019 / EN ISO 15223-
1:2016 / EN ISO 14971:2012 / ISO
13485:2016 / ISO 9001:2015

