

Protective Coverall Type 5-6



Product Code: **PPE002**

Description

40 gsm PP / SMS fabric - Sewn from light , soft , water resistant two layer non-woven fabric

Elasticated hood, wrist and ankle seams for maximum comfort

Can decrease the potential for cross - contamination in pandemic responseactivities when used with mask and gloves

Available in White colour

Certifications

ISO 13485:2016 / ISO 9001:2015(QMS)
/ CE ATTESTATION OF CONFORMITY
RELATED TO MEDICAL DEVICES
93/42/EEC/CLASS 1

