

# Protective Coverall Type 5-6



Product Code: **PPE001**

## Description

57 gsm NW textile backsheet white - Sewn from light, soft, water resistant two layer non-woven fabric

Elasticated hood, wrist and ankle seams for maximum comfort

Can decrease the potential for cross - contamination in pandemic response activities when used with mask and gloves

Available in White colour

## Certifications

EN 13795-1:2019 / EN ISO 15223-1:2016 / EN ISO 14971:2012 / ISO 13485:2016 / ISO 9001:2015

